

[View on rigaku.com](https://www.rigaku.com)

EDXRF1592 - Complex Lube Oils



Scope

The measurement of Mg, Si, P, S, Cl, Ca, Zn, Ba, and Mo in new lube oils is demonstrated using the empirical approach.

Background

Quality control and quality assurance during the lube oil blending and manufacturing process is essential. Basic lube oil formulations containing only P, S, Ca and Zn are often augmented with additional elements and additives to create lubricating oils for specific needs. A fast, simple method of analyzing lube oils is important throughout the QC/QA process. Rigaku meets this industry need with a high performance, low cost benchtop EDXRF system. Rugged and reliable, the [NEX QC+](#) is an ideal tool for elemental analysis, with simple and intuitive software for the non-technical operator.

Calibration

To demonstrate a complex fresh lube oil formulation, empirical calibrations were built using a suite of 23 commercially available calibration standards. The suite of calibration standards must be representative of the lube oil formulation to be analyzed. Use of empirical calibration ensures the particular lube oil formulation is exactly characterized and modeled, which yields optimum accuracy. Elements in the lube oil should span over each concentration range of interest, and the elements in the oil should vary independently of each other. Alpha corrections are then employed to automatically compensate for variations in X-ray absorption and enhancement effects within the sample due to the independent variations in element concentration. A summary of the empirical calibrations is shown here.

Element	Concentration range
Mg	0.020 – 0.080 %

Si	0.010 – 0.045 %
P	0.002 – 0.250 %
S	0.050 – 2.500 %
Cl	0.001 – 0.200 %
Ca	0.001 – 0.500 %
Ba	0.005 – 0.040 %
Zn	0.001 – 0.250 %
Mo	0.005 – 0.050 %

Repeatability

Representative samples from the calibration were chosen to demonstrate typical instrument repeatability (precision). Ten repeat analyses of each sample were performed with the sample in static position.

Sample: 7				
Units: %				
Element	Standard value	Average value	Std. dev	% Relative dev
Mg	0.025	0.0326	0.0017	6.8%
Si	0.030	0.0348	0.0016	5.3%
P	0.100	0.1007	0.0007	0.7%
S	1.251	1.272	0.007	0.6%
Cl	0	---	---	---
Ca	0.350	0.3578	0.0015	0.4%
Ba	0.010	0.0096	0.0003	3.0%
Zn	0.001	0.0011	0.0001	10%
Mo	0.015	0.0156	0.0001	0.7%

Sample: 9				
Units: %				
Element	Standard value	Average value	Std. dev	% Relative dev
Mg	0	---	---	---
Si	0.021	0.0201	0.0013	6.2%
P	0	---	---	---
S	0.750	0.7479	0.0034	0.5%

Cl	0.006	0.0053	0.0002	3.3%
Ca	0.300	0.3055	0.0017	0.6%
Ba	0.030	0.0305	0.0003	1.0%
Zn	0	---	---	---
Mo	0.050	0.0520	0.0005	1.0%

Sample: 21				
Units: %				
Element	Standard value	Average value	Std. dev	% Relative dev
Mg	0.060	0.0695	0.0035	5.8%
Si	0	---	---	---
P	0.250	0.2492	0.0041	1.6%
S	2.500	2.514	0.014	0.6%
Cl	0.120	0.1181	0.0013	1.1%
Ca	0.010	0.0120	0.0002	2.0%
Ba	0	---	---	---
Zn	0.050	0.0487	0.0002	0.4%
Mo	0	---	---	---

Sample: 6				
Units: %				
Element	Standard value	Average value	Std. dev	% Relative dev
Mg	0	---	---	---
Si	0	---	---	---
P	0.010	0.0110	0.0003	3.0%
S	0.250	0.246	0.001	0.4%
Cl	0.180	0.1834	0.0018	1.0%
Ca	0.001	0.0019	0.0001	10%
Ba	0.003	0.0031	0.0001	3.3%
Zn	0.225	0.2198	0.0006	0.3%
Mo	0.003	0.0030	0.0001	3.3%

Conclusion

The NEX QC+ offers analysts and technicians a simple yet powerful and versatile system for quantifying elemental composition of new, fresh lube oils using the empirical approach. The results of this study indicate that given stable samples, proper sample handling and proper calibration technique, the Rigaku NEX QC+ EDXRF can achieve excellent results measuring the elemental concentration of lubricating oils. The performance shown in this and other App Notes NEX QC and NEX QC+ can be used to comply with ASTM D6481 for the measurement of P, S, Ca and Zn.

Related products



NEX QC Series

Combines quality, affordability, and performance for a wide range of applications