

EDXRF1366 - Analysis of Metals in Nickel Ore



Scope

This application note uses [NEX QC+](#) to demonstrate the analysis of important metal oxides in nickel ore using the empirical approach.

Background

During site characterization and preparation for smelting, it is important to determine concentrations of various metals in nickel ores. While NiO and Fe₂O₃ are the most important metal oxides, TiO₂, Cr₂O₃, and MnO may also be present in significant quantities. Characterizing the ore is important to determine potential yield, as well as during extraction and smelting. EDXRF is an excellent tool for mine site exploration and characterization, and can be used throughout the processing steps, extraction, and smelting to measure the ore material, concentrate, matte, filter cakes, and slags.

Calibration

20 site-specific calibration standards were supplied as pressed pellets for empirical calibration. Optimum calibration can be achieved using standards that represent the ore material at the site, where each standard is assayed for all oxides of interest and concentrations evenly span each range and vary independently. Given assays for light element oxides such as Al₂O₃, SiO₂, and SO₃, those calibrations can be made as well.

Component	Concentration range (mass%)	Standard error of estimate
NiO	0.06 – 3.00	0.0146
Fe ₂ O ₃	8 – 66	0.2306

TiO ₂	0.01 – 2.41	0.0184
Cr ₂ O ₃	0.04 – 3.64	0.0296
MnO	0.10 – 2.00	0.0157

Standard sample recovery and precision

Representative low, medium, and high NiO concentration calibration standards were measured 10 consecutive times each in a static position to demonstrate effective recovery and analytical precision.

Sample ID: Sample ST14			
Units: Mass%			
Oxide	Assay value	NEX QC+ average value	Std. dev
NiO	2.98	3.00	0.003
Fe ₂ O ₃	21.71	21.89	0.02
TiO ₂	---	0.019	0.004
Cr ₂ O ₃	1.07	1.106	0.007
MnO	0.38	0.366	0.002

Sample ID: Sample S103			
Units: Mass%			
Oxide	Assay value	NEX QC+ average value	Std. dev
NiO	1.28	1.297	0.006
Fe ₂ O ₃	65.63	66.19	0.07
TiO ₂	0.05	0.051	0.003
Cr ₂ O ₃	3.64	3.636	0.015
MnO	1.16	1.178	0.008

Sample ID: Sample 1696			
Units: Mass%			
Oxide	Assay value	NEX QC+ average value	Std. dev
NiO	0.73	0.746	0.003
Fe ₂ O ₃	12.21	12.16	0.02
TiO ₂	0.18	0.166	0.004
Cr ₂ O ₃	0.32	0.322	0.003
MnO	0.20	0.190	0.001

Sample ID: Sample 1702			
Units: Mass%			
Oxide	Assay value	NEX QC+ average value	Std. dev
NiO	0.26	0.257	0.001
Fe ₂ O ₃	12.14	12.05	0.02
TiO ₂	0.30	0.292	0.005
Cr ₂ O ₃	0.23	0.238	0.002
MnO	0.20	0.204	0.001

Analysis of unknown samples

Unknown samples with expected values were submitted from the mine site for measurement against the calibrations built. Results are shown below.

Sample ID	NiO		Fe ₂ O ₃		TiO ₂		Cr ₂ O ₃		MnO	
	Expected Mass%	Measured Result	Expected Mass%	Measured Result	Expected Mass%	Measured Result	Expected Mass%	Measured Result	Expected Mass%	Measured Result
ST11	2.41	2.43	24.95	25.27	<0.01	0.03	1.21	1.31	0.38	0.38
ST16	2.04	2.05	15.37	15.32	<0.01	ND	1.23	1.22	0.22	0.25
ST20	1.79	1.76	65.45	64.96	<0.01	0.05	2.73	3.08	1.22	1.24
AP06	1.33	1.41	16.38	17.17	0.02	0.03	0.60	0.63	0.27	0.27

Conclusion

The NEX QC+ offers analysts a simple yet powerful and versatile system for quantifying elemental composition using the empirical approach. The results of this study indicate that given stable samples, proper sample handling and proper calibration technique, the Rigaku NEX QC+ EDXRF can achieve excellent results for the measurement of the key metal oxides NiO, Fe₂O₃, TiO₂, Cr₂O₃ and MnO in nickel ores.

Related products



NEX QC Series

Combines quality, affordability, and performance for a wide range of applications