

[View on rigaku.com](http://www.rigaku.com)


B-XRD1113 - Evaluation of oxidation state by the BVS method

Introduction

The valence of metal ions in a crystal is an important structural parameter for understanding a material's property and for designing new materials. There are numerous methods of experimentally determining the valence of a metal ion, including XAFS (X-ray Absorption Fine Structure), XPS (X-ray Photoelectron Spectroscopy), Mössbauer absorption spectroscopy, MEM (Maximum Entropy Method) and BVS (Bond Valence Sum). BVS is an empirical method based on Rietveld analysis that estimates the valence states based on the bond distances to neighboring atoms. Here, we illustrate the oxidation states of Mn oxides by BVS.

Measurements and results

X-ray powder diffraction of Mn_3O_4 and Mn_2O_3 were performed using Cu radiation. X-ray fluorescence of Mn is removed by the XRF reduction mode of the detector. The d-spacing is measured to be 0.82 \AA ($2\theta/\theta = 140^\circ$), which corresponds to half to a third of the length of the Mn-O bond distance. BVS, which is computed from the Mn-O bond distances obtained from the Rietveld analysis, and the Bond Valence Parameters ($\text{Mn}^{2+}\text{-O}^{2-}$: 1.790 (3) and $\text{Mn}^{3+}\text{-O}^{2-}$: 1.760 (5)) (1), predicts the proper oxidation state of the Mn ion for a given crystal structure as shown in Fig.1. Table 1 and 2 show the BVS values for Mn ions in the Mn_3O_4 and Mn_2O_3 crystal structures, respectively. The reliability factors of the Rietveld analysis are $R_{\text{wp}} = 4.49\%$, $R_p = 2.62\%$, $S = 1.6093$ for Mn_3O_4 , and $R_{\text{wp}} = 4.55\%$, $R_p = 3.59\%$, $S = 1.1949$ for Mn_2O_3 , respectively. BVS is useful for evaluating the stability and the validity of artificially designed crystal structures, in addition to being a simple approach for predicting and estimating the oxidation states of metal ions from the crystal structure.



The screenshot shows the 'Crystal Structure Information' window with the following input parameters:

- Space group: 141 : I41/amd:1
- a, Å: 5.76227; b, Å: 5.76227; c, Å: 9.47347
- α , °: 90.000; β , °: 90.000; γ , °: 90.000

The table below shows the data for the atoms/molecules:

Atom/Molecule	Density Peak	Distance	Angle	Element	x	y	z	Occ.	B	Valence	Bond Valen...
Mn1(Mn)				Mn	0.000000	0.000000	0.000000	1.000	0.13(3)	3	1.92276015...
Mn2(Mn)				Mn	0.000000	0.250000	0.625000	1.000	0.17(2)	3	2.97990813...
O1(O)				O	0.000000	0.2228(4)	0.3825(2)	1.000	0.19(4)	0	0

Figure 1: The input screen for the BVS calculation of Mn_3O_4 . (When the Mn valences are input as trivalent, the BVS of Mn1 and Mn2 sites are calculated as 1.92 (divalent) and 2.98 (trivalent), respectively.)

Table 1: BVS of Mn ions and structural parameters of Mn_3O_4

Space group		141: I4₁/amd:1			
a (Å)	5.76227(10)	b (Å)	5.76227(10)	c (Å)	9.47347(18)
α (°)	90.000	β (°)	90.000	γ (°)	90.000

Element	x	y	z	Occupancy	B	BVS
Mn1	0.0000	0.0000	0.0000	1.0	0.13(3)	2.085
Mn2	0.0000	0.2500	0.6250	1.0	0.17(2)	2.980
O1	0.0000	0.2228(4)	0.3825(2)	1.0	0.19(4)	–

Table 2: BVS of Mn ions and structural parameters of Mn₂O₃

Space group		206: Ia3			
a (Å)	9.41533(13)	b (Å)	9.41533(13)	c (Å)	9.41533(13)
α (°)	90.000	β (°)	90.000	γ (°)	90.000

Element	x	y	z	Occupancy	B	BVS
Mn1	0.0000	0.0000	0.0000	1.0	0.46(3)	3.047
Mn2	0.28519(6)	0.0000	0.0000	1.0	0.50(3)	3.026
O1	0.1296(3)	0.1473(3)	-0.0841(2)	1.0	0.50(3)	–

References

(1) I. D. Brown and D. Altermatt: *Acta Cryst.*, **B41** (1985) 244-247.

Related products



SmartLab

Advanced state-of-the-art high-resolution XRD system powered by Guidance expert system software



SmartLab SE

Highly versatile multipurpose X-ray diffractometer with built-in intelligent guidance



SmartLab Studio II

Windows-based software suite for Rigaku's X-ray diffractometers