

[View on rigaku.com](https://www.rigaku.com)

# B-XRD2027 - Evaluation of uniformity of a single crystal substrate by rocking curve measurement


## Introduction

Crystallographic non-uniformity of single crystal substrates and buffer layers, such as defects and orientation distribution, may decrease the performance of thin-film devices. XY mapping of the rocking curve allows us to evaluate the in-plane uniformity of single crystal substrates and buffer layers visually and quantitatively.


## Measurements and results

MgO 002 (diffraction angle  $2\theta=43.08^\circ$ ) rocking curve measurement ( $\omega$  scan) was performed on a 10 mm  $\times$  10 mm MgO substrate, with 0.5 mm intervals in the X and Y directions using a 0.5 mm  $\phi$  incident X-ray beam.

Figure 1 shows the color map of the position of the most intense peak observed in each profile. A change in peak position was observed over the entire MgO substrate. This indicates that the inclination of the MgO 002 lattice plane was non-uniform across the sample. Figure 2 compares profiles collected at the measurement positions marked with the dashed lines in Figure 1 (X = -2.0 mm, Y = -4.5 to +5.5 mm). The profiles at Y = +1 to +2.5 mm, which corresponds to the "red" part in Figure 1 (a), showed  $0.055^\circ$  higher-angle peak positions than other parts. This indicates that the MgO 002 lattice planes of this section tilt against the other parts of the sample.



**Figure 1:** Color maps of MgO 002 rocking curve peak position. (a) Full range map, (b) emphasizing medium range



**Figure 2:** MgO 002 rocking curve profiles collected at  $X = -2.0 \text{ mm}$ ,  $Y = -4.5$  to  $5.5 \text{ mm}$ .

---

## Related products



### SmartLab

Advanced state-of-the-art high-resolution XRD system powered by Guidance expert system software



### SmartLab SE

Highly versatile multipurpose X-ray diffractometer with built-in intelligent guidance



### SmartLab Studio II

Windows-based software suite for Rigaku's X-ray diffractometers