

[View on rigaku.com](https://www.rigaku.com)

# EDXRF1994 - Material Identification



## Scope

Demonstrate Material Identification for standard ASTM and JIS alloy tables and user-definable custom identification tables.

## Background

Material identification, also called positive material identification or PMI, is used to analyze and identify materials for quality assurance and quality control to verify proper alloy or material is being used. Material identification can also be used to quickly sort unknown alloys and other materials, including powders and liquids, and can be used in industrial forensics and failure material analysis identification. By looking at the structure of the elemental constituents in a sample, material identification verifies material type by comparing differences in composition that define the specific types of materials being analyzed.

XRF offers a non-destructive means for material identification. Applied Rigaku Technologies EDXRF systems offer the standard alloy identification tables as well as the capability for the user to create a custom table defined for any material where the types of materials can be differentiated by differences in element composition.

## Standard alloy identification tables

The ASTM and JIS standard alloy tables for ferrous, aluminum, and copper alloys are included.

Material Identification Result				
Standard		[All]		
Category		[All]		
No.	Matching	Standard	Category	Material type
1	98.10	JIS	for special purposes	SUS 316
2	98.10	JIS	for piping	SUS 316 TPY
3	98.10	ASTM	for special purposes	S 31600
4	98.10	ASTM	for special purposes	S 31603
5	98.10	ASTM	for special purposes	S 31635
6	98.10	ASTM	for piping	A 358 316
7	98.10	ASTM	for piping	A 358 316 L
8	98.10	ASTM	for piping	A 358 S 31635
9	98.10	ASTM	for piping	A 409 TP 316
10	98.10	ASTM	for piping	A 409 TP 316 L
11	98.09	JIS	for piping	SUS 316 TP
12	98.09	JIS	for piping	SUS 316 TiTP
13	98.09	JIS	for heat transfer	SUS 316 TB
14	98.09	JIS	for heat transfer	SUS 316 TiTB
15	98.09	JIS	for heat transfer	SUS 316 TF
16	98.09	JIS	for structural purposes	SUS 316 TKA
17	98.09	JIS	for structural purposes	SUS 316 TKC
18	98.09	JIS	Steel forgings	SUS 316

In this example, an unknown piece of stainless steel is measured using semi-quant FP and identified as SS316.

Material Identification Result				
Standard		[All]		
Category		[All]		
No.	Matching	Standard	Category	Material type
1	99.13	JIS	Copper alloy Castings	YB3C 3
2	98.77	JIS	Copper and copper al...	C 3601 B
3	98.77	JIS	Copper and copper al...	C 3602 B
4	98.77	JIS	Copper and copper al...	C 3601 W
5	98.77	JIS	Copper and copper al...	C 3602 W
6	98.67	JIS	Copper and copper al...	C 3603 B
7	98.67	JIS	Copper and copper al...	C 3604 B
8	98.67	JIS	Copper and copper al...	C 3603 W
9	98.67	JIS	Copper and copper al...	C 3604 W
10	98.67	ASTM	Copper and copper al...	C 37700
11	98.59	JIS	Copper and copper al...	C 3560 P
12	98.59	JIS	Copper and copper al...	C 3560 R
13	98.59	ASTM	Copper and copper al...	C 35600
14	98.49	JIS	Copper and copper al...	C 3771 B
15	98.48	JIS	Copper and copper al...	C 3561 P
16	98.48	JIS	Copper and copper al...	C 3561 R
17	98.46	ASTM	Copper and copper al...	C 36000
18	98.40	JIS	Copper and copper al...	C 3713 P

Identification of unknown brass sample

No.	Matching	Standard	Category	Material type
1	99.55	JIS	Aluminium and alumin...	A 6 N 01 S
2	99.41	JIS	Aluminium and alumin...	A 6151 FD
3	99.41	JIS	Aluminium and alumin...	A 6151 FH
4	99.41	ASTM	Aluminium and alumin...	6151
5	99.31	JIS	Aluminium and alumin...	A 3105 P
6	99.31	ASTM	Aluminium and alumin...	3105
7	99.17	JIS	Aluminium and alumin...	A 6063 BE
8	99.17	JIS	Aluminium and alumin...	A 6063 TE
9	99.17	JIS	Aluminium and alumin...	A 6063 TD
10	99.17	JIS	Aluminium and alumin...	A 6063 S
11	99.17	ASTM	Aluminium and alumin...	6063
12	99.14	JIS	Aluminium and alumin...	A 6061 P
13	99.14	JIS	Aluminium and alumin...	A 6061 BE
14	99.14	JIS	Aluminium and alumin...	A 6061 BD
15	99.14	JIS	Aluminium and alumin...	A 6061 W
16	99.14	JIS	Aluminium and alumin...	A 6061 TE
17	99.14	JIS	Aluminium and alumin...	A 6061 TD
18	99.14	JIS	Aluminium and alumin...	A 6061 S

Identification of unknown aluminum sample

## Custom identification table

Rigaku Material Identification software allows the user to create and define a custom identification table for any material where the types of the material are differentiated by elemental composition.

In this example, an identification for sorting unusual auto cats is demonstrated. The material types and criteria for differentiating are determined by the user to best fit the specific needs at the facility.

Create a new standard table, define the category, and list the major elements in the material that are used to differentiate the various material types.

**Standard Information**

Standard: Auto Cats Sort

Information: Unusual auto cats

OK Cancel

**Category Information**

Category: Auto Cats

Information:

Components:

Component	Unit	Summation
1 SiO2	mass%	
2 TiO2	mass%	
3 PbO	mass%	
4 ZrO2	mass%	
5 BaO	mass%	
6 La2O3	mass%	
7 CeO2	mass%	
8 MgO	mass%	
9 Al2O3	mass%	
10		

Insert Delete

OK Cancel

Next, list the types of material and enter the elemental limits that define the differences between each type.

Standard: Auto Cats Sort		Edit													
Category: Auto Cats		Edit													
Type	Note	SiO2 mass%		TiO2 mass%		PbO mass%		ZnO2 mass%		BaO mass%		La2O3 mass%			
		Lower limit	Upper limit	Lower limit	Upper limit	Lower limit	Upper limit	Lower limit	Upper limit	Lower limit	Upper limit	Lower limit	Upper limit	Lower limit	Upper limit
1 Diesel AITi		No limit		40	1	100	No limit	0.79	No limit	0.2	No limit	0.1	No limit	0.1	No limit
2 Diesel SiC		70	100	No limit		0.5	No limit	0.79	No limit	0.1	No limit	0.1	No limit	0.1	No limit
3 High Pb		No limit	40	No limit		0.5	0.8	100.0	No limit	20	0.15	2.00	No limit	1	No limit
4 Mixed Ceramics		No limit	40	No limit		0.5	No limit	0.79	No limit	20	0.15	2.00	No limit	1	No limit
5															
6															

Standard: Auto Cats Sort		Edit																			
Category: Auto Cats		Edit																			
Type	Note	PbO		ZnO2		BaO		La2O3		CaO2		MgO		Al2O3							
		Upper limit	Lower limit	Upper limit	Lower limit	Upper limit	Lower limit	Upper limit	Lower limit	Upper limit	Lower limit	Upper limit	Lower limit	Upper limit	Lower limit						
1 Diesel AITi		0.79	No limit	0.2	No limit	0.1	No limit	0.1	No limit	0.1	No limit	0.1	No limit	10	No limit	40					
2 Diesel SiC		0.79	No limit	0.1	No limit	0.1	No limit	0.1	No limit	10	No limit	10	No limit	40	No limit						
3 High Pb		100.0	No limit	20	0.15	2.00	No limit	1	1	10	No limit	10	No limit	40	No limit						
4 Mixed Ceramics		0.79	No limit	20	0.15	2.00	No limit	1	1	10	No limit	10	No limit	40	No limit						
5																					

Material Identification Result

Standard: [All]

Print

Category: [All]

No.	Matching	Standard	Category	Material type
1	99.25	Auto Cats Sort	Auto Cats	Diesel SiC

Example of identification of SiC diesel auto cat sample

Material Identification Result

Standard: [All]

Print

Category: [All]

No.	Matching	Standard	Category	Material type
1	98.03	Auto Cats Sort	Auto Cats	Diesel AITi
2	86.96	Auto Cats Sort	Auto Cats	Mixed Ceramics
3	86.17	Auto Cats Sort	Auto Cats	High Pb

Identification of an aluminum titanate diesel cat.

No.	Matching	Standard	Category	Material type
1	98.05	Auto Cats Sort	Auto Cats	High Pb
2	97.68	Auto Cats Sort	Auto Cats	Mixed Ceramics
3	95.66	Auto Cats Sort	Auto Cats	Diesel AITi

Identification of an auto cat with unusually high lead

## Conclusion

Rigaku Material Identification software gives the user a tool for alloy PMI and the ability to create custom identification tables for material verification and sorting applications.

---

## Related products



### NEX CG II Series

High-performance *indirect excitation* EDXRF for complex applications with trace elements and variable base matrices



### NEX DE Series

High-power 60 kV EDXRF systems delivering speed, precision, and small spot measurements



### NEX QC Series

Combines quality, affordability, and performance for a wide range of applications