

B-XRD1089 - Micro-area mapping measurement of printed circuit boards

Introduction

When analyzing a sample consisting of an extremely large number of micro-size parts, such as PCBs (printed circuit boards), it is necessary for the evaluation of the crystalline state of each part to use a focused X-ray beam and irradiate only a specific part. Optical systems in the past employed slits and collimators to convert the divergent beam emitted from the X-ray focal spot to a point focus beam. The downside to this technique is the large attenuation of the X-ray beam intensity, requiring a long time for micro-area measurement to be completed. By using the SmartLab μ HR, a SmartLab model with microfocus optics, in combination with a high-speed 2D X-ray detector, a high-intensity micro-area measurement is possible in a short time even though the XG output is 1/10 of those in the past.

Measurements and results

The SmartLab μ HR is equipped with a microfocus rotating anode X-ray tube. The X-ray beam is monochromatized and converted to a parallel beam 1 mm $\times$ 1 mm in diameter by using CMF microfocus optics. The beam width is narrowed by a collimator. As a result, as indicated in Fig. 1 and Table 1, compared with a collimator and the CBO-f point focus optics used in the past, a high-resolution and high-brilliance measurement was achieved even though the XG output is less than 1/10.

Micro-area measurement of PCB was performed using the SmartLab μ HR. In both measurement points, identification of compounds included in micro-area of PCB and the orientation of coarse particles can be confirmed in an extremely short amount of time of approx. two minutes.

Table 1: Intensity and FWHM of optics

Optics	XG output (kW)	Beam size (mm)	Normalized intensity	Brilliance (Intensity / Beam size)
SmartLab μ HR	0.8	0.3 ϕ	52.10	579
SmartLab CBO-f	9.0	0.5 ϕ	53.45	68.1
SmartLab collimator	9.0	0.3 ϕ	1.00	3.5




Figure 1: Data comparison of optics (Sample: Si)




Figure 2: Micro-area measurement of PCB (Measurement time approx. two minutes)

Related products



HyPix-3000

Compact photon counting x-ray detector



SmartLab

Advanced state-of-the-art high-resolution XRD system powered by Guidance expert system software



SmartLab Studio II

Windows-based software suite for Rigaku's X-ray diffractometers