

PX026 - Sulfur SAD phasing of thaumatin in 10 minutes with the HyPix-Arc 150° and XtaLAB Synergy-DW

Introduction

T. danielli thaumatin, a sweet-tasting protein often used as a sweetener, features an axis of ~ 150 Å. It is well characterised by protein crystallography and makes a good candidate for testing new equipment as it becomes available. Recently Rigaku Oxford Diffraction introduced a new detector: the HyPix-Arc 150°. The new detector provides extremely high theta coverage which can be used for fast collection of Friedel pairs. The detector features HPC technology for the highest accuracy of intensity measurement for both weak and strong signal ensuring that anomalous signal is well measured for S-SAD phasing applications. This note demonstrates fast S-SAD phasing of *T. danielli* thaumatin.

Experimental

T. daniellii thaumatin was purchased from Sigma and suspended in ultrapure, laboratory grade water to a concentration of 50 mg/ml. Crystals of thaumatin were grown at 18°C from hanging drops consisting of 6 µL of 50 mg/ml thaumatin plus 4 µL of crystallization solution [24% (w/v) NaK Tartrate, 15% (v/v) ethylene glycol, and 0.1 M BisTris Propane, pH 6.6] over a well containing 500 µL of crystallization solution. Prior to flash-cooling, each crystal was soaked for ~ 15 to 30 seconds in the above crystallization solution at 25% (v/v) ethylene glycol for cryoprotection. A crystal was then selected and stored in a dry shipper prior to data collection.

Table 1: XtaLAB Synergy-DW with HyPix-Arc 150° key features

X-ray source	Type	PhotonJet-DW dual wavelength rotating anode
	Power	1.2 kW
	Divergence control	Continuously variable slit
Goniometer	Type	4-circle Kappa goniometer
	XTD range	30-250 mm
Detector	Type	HyPix-Arc 150° hybrid photon counting detector
	Pixel size	100 µm




Figure 1: Thaumatin crystal used.

The data were collected on a XtaLAB Synergy-DW dual wavelength rotating anode diffractometer equipped with the curved HyPix-Arc 150° detector. A crystal measuring 0.13 x 0.15 x 0.27 mm was selected and screened with CrysAlisPro to determine it had a suitable diffraction limit. Using CrysAlis^{Pro} a Cu strategy was planned that aimed to get suitable data for structure solution using S-SAD phasing in under 10 minutes. Due to the superior spatial resolution of the 100 μm pixels with no point spread, it was possible to collect the data at a crystal to detector distance of 60 mm while still resolving the long axis. The strategy module was also set to align the crystal for best orientation of the long axis on the detector. The data that were collected were processed in CrysAlis^{Pro} with a 1.45 Å resolution limit and passed on to AIMLESS for merging. The Phaser-EP pipeline in the CCP4 suite was used to solve the structure. PARROT and BUCCANEER were then used for density modification and model building respectively and Coot was used for model correction.

Results

The sample diffracted extremely well and for the chosen scan speed of 1.5 seconds per degree (0.3 s per 0.2° image), a diffraction limit of 1.45 Å was determined. The unit cell was found by CrysAlis^{Pro} to be tetragonal P with $a = b = 57.81$ Å, $c = 149.95$ Å and the cell was found automatically with no user intervention required.

Divergence slits were fully open throughout the experiment and despite this, the long axis was well resolved at an XTD of 60 mm as can be seen in figure 2.

Table 2: Dataset statistics

Space group	P4 ₁ 2 ₁ 2
Unit cell (Å)	57.81, 57.81, 149.95
Total time	9m 53s
Resolution (Å) [last shell]	22.96 - 1.45 [1.48 - 1.45]
Total / unique reflections	699245 / 45973
Completeness [last shell]	99.7 [94.3]

Multiplicity [last shell]	15.2 [3.3]
$\langle I/\sigma(I) \rangle$ [last shell]	23.1 [1.4]
R_{merge} (%) [last shell]	7.5 [69.6]
$CC_{1/2}$ [last shell]	0.999 [0.571]
$R1 / R_{\text{free}}$ (%)	12.0 / 16.0




Figure 2: Peaks along the long 150 Å axis are well resolved at 60 mm

This excellent spatial resolution can be attributed to the 100 µm pixel size and the single-pixel top-hat point spread function.

After processing with CrysAlis^{Pro} and merging with AIMLESS an R_{merge} of 7.5% was obtained to 1.45 Å. Solution with the Phaser pipeline in CCP4 gave 24 initial sites from which the full model was easily built with automated building routines in Buccaneer. Following iterative cycles of refinement, and model correction with COOT, the final model produced $R1 = 12.0\%$ and $R_{\text{free}} = 16.0\%$.





Figure 3: Nicely defined electron density map following S-SAD phasing from a 10 minute thaumatin dataset

Conclusion

The crystal structure of *T. daniellii* thaumatin was easily solved by S-SAD phasing to 1.45 Å. The extremely high data quality provided by the HyPix-Arc 150° and XtaLAB Synergy-DW diffractometer allowed this to be achieved with only 10 minutes of data collection time.

Related products



HyPix-Arc 150°

A curved detector based on direct X-ray detection technology with the highest 2θ range at a single position available for the home lab.



XtaLAB Synergy-DW

Spectacular performance combined with dual wavelength versatility, provides the perfect answer for high-capacity Chemical Crystallography labs or for X-ray facilities that support Chemical Crystallography and Protein Crystallography.