

# B-TA1043 - Evaluation of CO<sub>2</sub> adsorption below ambient temperature by low-temperature STA

## Introduction

Gas adsorption on a substance varies with temperature. Since the amount of adsorption generally increases at lower temperatures, it is desirable to evaluate the amount of adsorption up to a temperature range below room temperature. In this study, low-temperature STA was used to measure the amount of CO<sub>2</sub> adsorbed on HKUST-1, MOF (metal-organic structures), below room temperatures to evaluate CO<sub>2</sub> adsorption versus temperature.

## Measurement and analysis example

The results of the low-temperature STA for sample HKUST-1 are shown in Figure 1. The STA equipped with a refrigerated cooling unit allows cooling down to -40°C and heating program measurements. The sample was pretreated by heating to 180°C in helium atmosphere, and then was cooled down to -40°C. The atmosphere was then switched to CO<sub>2</sub> at each target temperature, and the amount of CO<sub>2</sub> gas adsorbed by the samples was measured. The CO<sub>2</sub> adsorption temperatures were kept at -29°C, -4°C, 12°C, and 23°C. After measuring the amount of adsorption at each temperature, the sample was reheated up to 180°C in helium atmosphere each time.





Figure 1: Results of low temperature STA measurement of HKUST-1



**Figure 2:** Plot of temperature vs. adsorption amount

The results of STA show a 40% mass loss during the first heating process. After cooling down and changing the atmosphere to CO<sub>2</sub> at -29°C, a 62% mass increase due to CO<sub>2</sub> gas adsorption is observed. Similarly, at -4°C, 42%, at 12°C, 31%, and at 23°C, 23% of CO<sub>2</sub> gas adsorption is observed, respectively. Figure 2 shows a plot of adsorption versus temperature obtained from a series of measurements. Thus, the low-temperature STA makes it possible to evaluate the temperature dependence of the adsorption properties of materials in the temperature range below room temperature.

### Recommended equipment and software

- ▶ [STA8122](#) + refrigerated cooling unit
- ▶ [Vullios](#) measurement and analysis software

(Samples provided by Prof R. Matsuda and Prof S. Kusaka, Nagoya University)