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
Rigaku
Rigaku Battery
Application Ebook

About Rigaku

Who We Are

Since its inception in 1951, Rigaku has been at the forefront of analytical and industrial instrumentation technology. Today, with hundreds of major innovations to their credit, the Rigaku group of companies are world leaders in the fields of general X-ray diffraction, thin film analysis, X-ray fluorescence spectrometry, small angle X-ray scattering, protein and small molecule

X-ray crystallography, Raman spectroscopy, X-ray optics, semiconductor metrology, X-ray sources, computed tomography, nondestructive testing and thermal analysis.

Corporate Mission

To contribute to the enhancement of humanity through scientific and technological development.

RIGAKU BATTERY APPLICATION EBOOK

Fuel Cell

Rigaku's analysis technologies support R&D for safer and more efficient cells.

Fuel Cell

Separator
Gas diffusion layer
Electrolyte membrane

Analyze: Electrolyte membrane
For: Evaluation of thermal stability
Thermal analysis

Catalyst

Analyze: Catalyst
For: Optimizing electrochemical performance and quality assurance
Particle size and structure analysis
Elemental analysis

MEA, Membrane Electrode Assembly (Gas diffusion layer, catalyst layer and electrolyte membrane)

Analyze: MEA, Membrane Electrode Assembly
For: Optimizing electrochemical performance
Non-destructive imaging

RIGAKU BATTERY APPLICATION EBOOK

Analysis of Adhesion within Both Sides of a Fuel Cell (PEFC) Electrode

Analyze: Fuel cell
Use: Quality assurance
Elemental Analysis

Analyzed materials:
Pt as catalyst

In a polymer electrolyte fuel cell (PEFC: Polymer Electrolyte Fuel Cell), a cathode layer and an anode layer are placed on both sides of the electrolyte membrane. Since the amount of platinum (Pt) catalyst on each electrode layer is different, each layer must be managed separately. Conventional control methods require each surface to be measured separately, resulting in a large workload and time cost. X-ray fluorescence analysis is an effective method for efficiently solving these problems because it has high penetrating power, is nondestructive, and can simultaneously evaluate the amount of Pt adhered on both surfaces in a single measurement.


Figure 1: Schematic diagram of measurement

Measurement surface	A-side	B-side
Analysis value	0.001	0.42
Design value	0.003	0.405
Relative error (%)	1.0	3.7

Table 1: Analysis results for each aspect

Conclusion: Both sides (A-side and B-side) of the electrolyte membrane sheet coated with platinum-loaded carbon electrode were measured for 60 seconds using X-ray fluorescence of platinum and quantitatively analyzed by the thin film FP method. As shown in Table 1, the analysis results for each measurement point agreed with the design values for each surface within 15% relative error, confirming that the platinum loading in this film formation process was properly controlled.

RIGAKU BATTERY APPLICATION EBOOK

Evaluation of Water Absorption and Swelling in Polymer Electrolyte Membranes

Analyze: Fuel cell
Use: Evaluation of thermal stability
TA: Thermal Analysis

Analyzed materials:
Electrolyte membrane

In the fuel cell operating environment, changes in temperature and humidity affect the expansion and contraction of materials. Dimensional changes due to water absorption in polymer electrolyte membranes can have a significant impact on cell sealing and adhesive durability. There are limited methods to evaluate these properties in a humid environment, and STA and TMA with controlled humidity are effective approaches in material design and selection.


Figure 1: (a) weight change and (b) dimensional change of polymer electrolyte membranes as a function of relative humidity.

Conclusion: STA measurements under humidity-controlled conditions allow quantitative evaluation of water absorption by measuring the weight change of a sample under each relative humidity condition. TMA can also be used to determine the swelling behavior under each humidity condition.

RIGAKU BATTERY APPLICATION EBOOK

Internal Structure Evaluation of Fuel Cell Membrane Electrode Assembly (MEA)

Analyze: Fuel cell
Use: Optimizing electrochemical performance
Non-destructive analysis

Analyzed materials:
MEA, Membrane Electrode Assembly

For improving the performance and quality of fuel cells, it is important to evaluate the internal structure of the gas diffusion layer (GDL) in the membrane electrode assembly (MEA), a multifunctional base material responsible for gas diffusion, electron collection, and discharge of produced water. GDLs are composed of low-density materials such as carbon fiber, and conventional X-ray CT imaging using high-energy X-rays does not provide contrast because X-rays are not absorbed by GDLs, making it difficult to analyze pore size and diffusion pathways. The nanoCT, which is capable of low-energy X-ray CT imaging, makes it possible to visualize and evaluate the structure of GDLs, which is expected.

Figure 1: X-ray CT image of a membrane electrode assembly (MEA) and analysis of gas diffusion layer (GDL) pores

Conclusion: CT imaging using low-energy characteristics: X-rays of the NanoCT enabled high-contrast and clear observation of GDLs even in the presence of electrolyte membranes with high-density platinum catalyst attached. By visualizing the internal structures of all the substrates comprising the prototype MEA in three dimensions and quantitatively evaluating structural indices such as the GDL pore size, it will be possible to study the optimization of the MEA preparation conditions.

RIGAKU BATTERY APPLICATION EBOOK

SAXS-RMC Method for Estimating 3D Structure of Pt/GDC Supported Nanoparticles

Analyze: Fuel cell
Use: Optimizing electrochemical performance
Particle size and structure analysis

Analyzed materials:
Pt as catalyst

Catalysts consisting of platinum nanoparticles supported on gadolinium-doped ceria (Gd-doped ceria: GDC) (Pt/GDC catalysts) have been reported to exhibit better four electron oxygen reduction reaction activity than commercially available Pt/C catalysts. It is important to ensure that the reactions at the catalyst surface in the catalyst bed are facilitated, e.g., by the oxygen and hydrogen gas supply pathways. To understand these properties, it is useful to evaluate the three-dimensional structure of the entire catalyst layer, which is a complex aggregation of nanoparticles. However, with conventional observation methods, it is difficult to nondestructively determine how the nanoparticles are aggregated in three dimensions. By combining small angle X-ray scattering (SAXS) and reverse Monte Carlo (RMC) methods, it is possible to estimate the three-dimensional aggregation structure of nanoparticles non-destructively and statistically.


Figure 1: The resultant structural models obtained by RMC to match the observed SAXS patterns. Red and gray spheres for Pt/GDC particles, respectively.

Figure 2: A typical TEM image of Pt/GDC.

Conclusion: The structural model constructed by the combination of SAXS and RMC methods was in good agreement with the TEM observation results. Furthermore, the pore size distribution calculated based on the model was in good agreement with the experimental results obtained by nitrogen gas adsorption, confirming the reproducibility of macroscopic structural indices. These results indicate that the SAXS-RMC method is a very useful 3D structural analysis technique that is superior in both visualization of local structure and statistical structure reproduction.


Figure 3: The SAH analysis result curves are shown in the dashed line with the measurement data points in open circles. Average curves of three RMC simulations run for SAXS patterns are shown in the solid lines.

© 2025 RIGAKU HOLDINGS CORPORATION AND ITS GLOBAL SUBSIDIARIES. ALL RIGHTS RESERVED.

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